

Features of Versa-Pac® Accessory Products



MV75210

Compression and Support

Versa-Lastic® bandages are a comfortable and convenient way of securing a gel pack to the affected area. The bandage provides compression while allowing freedom of movement.

Compression from the Versa-Lastic® bandage provides even pressure around the injured area and offers support. A healing muscle injury can benefit from the use of a protective and supportive elastic bandage which may aid in the prevention of a re-injury.

Versa-Lastic® elastic bandages are constructed of a durable **LATEX-FREE** stretchable white fabric offering excellent memory and reusability. Easy to use sewn-in velcro closures eliminate the need for hooks or pins. Available in a variety of styles to address many treatment areas and options.

REF MV75100	3 in. x 48 in.	1/pkg
REF MV75110	6 in. x 48 in.	1/pkg
REF MV75150	3 in. x 24 in.	2/pkg

REF MV75140
This kit contains: 1 set/pkg
1-6 in. x 48 in. wrap with built-in pocket and
1-5 in. x 10.5 in. reusable hot/cold gel pack

Protect The Skin

To prevent the possibility of burns or frostbite, hot and cold packs should never be placed directly on the skin. A protective insulator, such as the Versa-Pac® reusable cloth cover should be used. This envelope style cover is custom sized to fit a variety of gel packs up to 6 in. x 11 in. The soft terry cloth cover provides patient comfort, insulation from the surface of the pack, and rapid temperature transfer.

REF MV75200 7 in. x 11.5 in. 6 pkg/case

REF MV75210
This kit contains: 1 set/pkg
1-7 in. x 11.5 in. cloth cover and
1-5 in. x 10.5 in. reusable hot/cold gel pack

Moist Therapy

Moistening the Versa-Pac® cloth cover accelerates temperature transfer and penetrates the injured tissue more quickly.

1313 West Grant Blvd. Fax 651-565-3971
Wabasha, MN 55981 Phone 651-565-2601
Toll free 1-800-328-9454 www.uni-patch.com
Email sales@uni-patch.com

© 2006 Tyco/Healthcare/Uni-Patch